Hera

HERAVEST M PRINT +

Ⓧ DE	Gebrauchsanweisung	3
Ⓧ GB	Instructions for use.....	5
Ⓧ FR	Mode d'emploi.....	7
Ⓧ IT	Istruzioni per l'uso	9

KULZER
MITSUI CHEMICALS GROUP

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Stand: 2021-10

Heravest M print+ is a phosphate bonded, graphite free and very fine grain precision investment compound especially for partial denture frameworks made with 'Rapid-Prototype' technologies. The investment material was specially developed for printed (Rapid Prototyping) and milled wax and acrylic frameworks. It is suitable for speed preheating and enables casting of all types of non-precious metal partial denture alloys (except titanium).

Storage

Keep the powder in a cool, dry place.
Do not store the mixing liquid at less than +5°C (41°F) as it is sensitive to low temperatures.

Working temperature

approx. 23°C (73°F) (Room temperature).

Mixing ratio

100g Powder – 20ml Liquid.

Expansion control

75% – 85% Liquid concentration
recommended Liquid concentration 80%

Concentrations			
Mixing ratio Powder/Liquid	Concentration (%)	Sachet 1 x 400g	
		Liquid (ml)	Distilled water (ml)
400g–80ml	75	60	20
400g–80ml	80	64	16
400g–80ml	85	68	12

The concentration may be adjusted as necessary.

Stirring

Use a spatula to stir the investment material by hand until the powder has been wetted thoroughly.

Mixing in vacuum

Place in vacuum for 15 sec but without switching the paddles on.
Mix for 60 sec. Check the vacuum/paddles frequently for correct functioning. Inadequate vacuum, results in an inaccurate fit and bubbles on the casting.

Working time

Approx. 5min at approx. 23°C (73°F) (Room temperature).
The working time depends on the room temperature. Warmth curtails the working time.

► Suitable Muffle Systems	
Ring free	Crucible former, die, size 1 and size. 2, Silicone sleeve, die, size 1 and size 2

Investing

Place the cast objects in the muffle mould and ensure that they are horizontal and slightly above the half muffle height. OK frameworks can be tipped slightly dorsally so that the air under the curved palatal plate can escape more easily during investment. The distance to the muffle wall and the upper surface of the muffle should be min. 10 mm.

Measure, mix and let Heravest M print + flow to the object without vibration as described above. Then switch on the vibrator with a high vibration frequency (6000 min⁻¹) but low intensity until the investment material is completely covered. Switch off the vibrator, fill the ring and do not vibrate again. **Important: Do not perform pressure embedding.**

Placing the rings

Roughen the top of the mould before placing it in the furnace!

Preheating Speed preheating

Precisely 20 min. after starting to mix the investment material, put the muffle with the funnel facing downwards in the oven preheated to 900 °C – max. 950 °C (1652 °F – max. 1742 °F) (please observe the precisely defined time, measured from the first contact between powder and liquid).

► Speed preheating			
Initial temperature	Final temperature	Retention time (min.) (Starting once the final temperature is reached again)	
		size 1	size 2
900°C - max. 950°C * (1652°F - max. 1742°F *)	900°C - 1030°C ** (1652°F - 1866°F **)	60	60

* If the required preheating temperature is above the max. initial temperature of 950 °C (1742°F), the muffle is put in the oven at 950 °C (1742°F) max. and heated to the required final temperature after 15 minutes.

** Observe the preheating temperature of the alloy. Note: Kulzer recommends in combination with Heravest M print + a preheating temperature of 950 °C (1742°F). That temperature is completely sufficient for mould filling as well as for the completed thermal expansion.

⚠ Safety Instructions

- Do not open the furnace during the wax burnout phase.
- Investments contain quartz. Do not inhale the dust! Risk of lung damage (silicosis, lung cancer). Recommendation: Wear a protective face mask, type FFP 2 EN 149. Prevent the release of dust when opening the bag and filling the mixing bowl.
- Rinse the empty bag with water before crumpling it up.
- Always moisten any dust before removing it from the workplace.
- After casting, place the completely cooled mould in water until it is fully moistened to prevent dust when devesting.
- Use an extractor with a fine dust filter when sandblasting.

We always attempt to include the latest results for our investment materials (development and application) in the instructions. Excellent precision of fit and surfaces are achieved with Kulzer equipment and materials. The results may vary if other materials (e.g. casting ring liners, muffle systems, acrylics, equipment etc.) or different techniques are used.

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Kulzer GmbH