

SEPTEMBER / OCTOBER HIGHLIGHTS 2024

durodent dental



NOW AVAILABLE THROUGH DURODENT
DENTSPLY DIGITAL RANGE

CAD CAM PMMA Milling Discs Denture Base and Denture Teeth / Crowns

**AVAILABLE
NOW**

**PURCHASE ANY 3 DISCS
INTRO SAVE 10%**

PLEASE NOTE:
SOME STOCK MAY REQUIRE LEADTIME
AND BE BACKORDERED INITIALLY

Multilayer PMMA - 20mm x 98mm

The Dentsply Sirona Multilayer PMMA Disc is now a validated dual-use material for temporary crowns and bridges and premium denture teeth.

- Enables labs to mill tooth arches or segments for complex denture cases
- Mills efficiently and polishes easily
- Composed of 13 shade layers that blend seamlessly together
- Provides natural fluorescence
- High flexural strength > 80 MPa
- All 16 VITA1 shades and one bleach shade are available
- Application: Premium full arch and tooth segments, Provisional Crowns and Bridges

**MULTILAYER
\$107.90 + GST**




Shaded (Mono) PMMA - 20mm x 98mm

Primarily for provisional use in complex restorations or when a long healing period is required before the permanent restoration is placed.

- Mills efficiently and polishes easily
- All 16 VITA1 shades and bleach are available
- Provides natural fluorescence
- High flexural strength > 80 MPa
- Fabrication with: inLab MC X5 and all open CAM-systems
- Application: Provisional Crowns and Bridges

**SHADED (MONO)
\$65.87 + GST**



Denture Base PMMA LUCITONE - 30mm x 98mm

Lucitone® is the industry benchmark for high strength, high-impact resistant denture base acrylics with over 40 years of clinical success.

- For full dentures including implant supported
- High impact and flexural strength, combined with excellent aesthetics, support denture production using the digital workflow
- Produce digital dentures recognized for precision & better fit
- Maintain high stability with 30 hours bur life
- Available in seven shades with and without veins



**LUCITONE DENTURE BASE
\$101.48 + GST**

Subscribe @ Durodent.com.au for our monthly email broadcast
Keep informed – New products, highlights, special offers

SUBSCRIBE

Valid SEPT / OCT 2024 - All Prices Exclude GST – Prices Do Not Include Freight - Durodent Reserves The Right To Amend Errors and Prices – Subject to stock availability

ONLINE STORE
www.durodent.com.au

SAVE TIME - ORDER ONLINE



Speak to our team:
Call us 1800 653 439

DIGITAL PRINT RESIN RANGE



LUCITONE 3D PRINT

**Dentsply Sirona Lucitone Digital Print™
3D Denture Base Resin**

Original Pink 906200 or Light Pink 906202

1 Bottle (1kg)

Validated printers for use with the Lucitone Digital Print Denture System: Carbon® M-Series, Asiga MAX™ UV, Asiga PRO 4K™. [1,2]

Validated post-processing units:
Dentsply Sirona inLab Speedcure® Processing Unit (REF #9494800 and 9494801)



\$669.50 + GST



**AVAILABLE
NOW**



dima®

1000g Bottles

Denture Teeth \$695.00 + GST

Denture Base \$650.00 + GST

Try In White \$540.00 + GST

Temp \$698.00 + GST

Splint Hard \$490.00 + GST



**BUY 2 + BOTTLES
SAVE 10%**

DETAX FREEPRINT®

500g and 1000g Bottles

CROWN Tooth Shade
From \$899.00 + GST

Denture Base
From \$306.00 + GST

Temp From \$365.00 + GST

Splint HARD
From \$288.00 + GST

Splint FLEX \$695.00 + GST



**BUY 2 + BOTTLES
SAVE 10%**



Economical Value Model Resins



PRECISE MODEL (ALMOND)

Easy-to-use, high-quality, and accurate in many printers

1L or 5L Bottle

MSLA and DLP Formula

**1L BOTTLE
\$145.00 + GST**



THERMO MODEL (GREY)

Higher Heat Resistance For Thermofforming Models

1L or 5L Bottle

MSLA and DLP Printer Formula

**5L ECO BOTTLE
\$435.00 + GST**

**BUY ANY 2 + BOTTLES
SAVE 10%**

Valid SEPT / OCT 2024 - All Prices Exclude GST - Prices Do Not Include Freight - Durodent Reserves The Right To Amend Errors and Prices - Subject to stock availability

ONLINE STORE
www.durodent.com.au

SAVE TIME - ORDER ONLINE



Speak to our team:
Call us 1800 653 439

DuraFlex Flexible Partial Denture Material – Reliance Dental USA

10% OFF

A partial denture material with flex that is suitable for short / medium term applications.

Pink Shade 99

(Not recommended for long term use in all applications)

Powder / Liquid Mix

Treated like a selfcure acrylic. See intructions at Durodent.com.au



Zhermack NEOCOLLIOD Alginate 500g C302205

High Precision Alginate

Normal Setting

High Flowability

3:30 (min:s)

\$21.00 ea + GST

**BUY 10 BAGS
GET 2 FREE**



Laboratory Kit
160g Pwd + 120ml Liq
\$178.20 ea + GST

Economy Kit
320g Pwd + 240ml Liq
\$315.00 ea + GST

Durodent Pros Stone (Yellowstone)

20kg

The Perfect Lab Stone

PSI 8,500 (58Mpa)
Set Time: 11–14 Mins
Exp: 0.12%
Ratio:30ml/100g



\$84.00 ea + GST

**BUY 4 BOXES
SAVE 10%**

Durodent Acrylpol 9kg Box

10% OFF

Pumice Like Acrylic Finishing Media - Pumice Alternative

Unique abrasives to cut faster and smoother than pumice. One grit size does it all.

No scratching or grooving.

L asts longer than pumice type products. No odor or bacteria. White.



\$103.50 + GST

ONLINE SHOP CLEARANCE OFFERS

OVERSTOCKED
END OF LINE
SHORT CODE
DISCONTINUED

www.durodent.com.au/clearance



Valid SEPT / OCT 2024 - All Prices Exclude GST – Prices Do Not Include Freight - Durodent Reserves The Right To Amend Errors and Prices – Subject to stock availability

ONLINE STORE
www.durodent.com.au

SAVE TIME - ORDER ONLINE



Speak to our team:
Call us 1800 653 439

LAB EQUIPMENT

Wassermann Wasi Steam CLASSIC II 174996A Manual Fill Steam Cleaner

LONG TIME TOP SELLER

ROBUST

STAINLESS STEEL TANK

METAL CONSTRUCTION

1 - 5 BAR ADJ PRESSURE

3.3L BOILER



\$2895.00 + GST

Wassermann Polishing Compact Unit WP-Ex 2000 II 146978

*Premium Quality
Engineered Product*


*2 Speed Motor - 1500 /
3000 RPM*

Integrated Daylight LEDs

680W x 486H x 527Dmm

*Suction with Dual Dust
Bags*

Made In Germany



\$7895.00 + GST
(RRP \$8968.00 + GST)

OMEC ASP/ PL4 Polishing Unit With Suction

10% OFF

Economical and Reliable

2 Speed Lathe Motor
1400-2800 RPM - 300W

Two Rubber Splash Backs
2 x Spindle Brush Chucks
36W Lighting System
Protective Screens
Single Suction Filter Bag

Made In Italy

510H x 670W x 500D

(Previously known as
Gamberini)



\$5355.00 + GST

Renfert PILLO 50224000 Pneumatic Deflasking Chisel

10% OFF

*Includes: Chisel Insert No. - Flat Chisel Wide 5022-0200
Air Consumption 9L/Min
Noise Level 79.3dB - Piston Freq 110Hz
Nominal Pressure 5 bar (72psi)
Working Pressure 4-6 bar
Weight 300g
3 Year Warranty*



\$595.00 + GST

Ultrasonic Cleaner CD4800

10% OFF

Digital Control

Plastic Basket

Cleaning Time:
5 Working Cycles

Ultrasonic Freq 42,000Hz

Tank Capacity 1.4L

160H x 180W x 230D
mm

Weight 1.75kg



\$207.00 + GST

Denstar Wax Induction Heater and Pot

Made In South Korea

10% OFF



**DS160N Wax
Instrument Heater
\$346.50 + GST**

**DS420 Wax Pot
1 Large Well
\$220.50 + GST**



Valid SEPT / OCT 2024 - All Prices Exclude GST - Prices Do Not Include Freight - Durodent Reserves The Right To Amend Errors and Prices - Subject to stock availability

SR Ivolen Tray - Ivoclar

A stable non-stick multi-purpose tray material characterised by its optimum setting phase.

It can be used with all dental impression materials applied with adhesives.

Colour: Yellow



1kg (2 x 500g)
Powder Only
\$130.50 + GST

1kg (2 x 500g) +
500ml Pack
\$180.00 + GST

Pala Putty – Kulzer

90 Shore A Hardness A
1:1 Mix

High-precision laboratory
A-silicone with a fast
setting time.

Resistant To Heat



3kg Pack
\$149.85 + GST

10kg Pack
\$396.00 + GST

Al Dente Ausblockwachs Thermoforming Blockout Wax - White 110g

581-004-011

For blocking out
undercuts prior to
Thermoforming



\$51.30 + GST

Kerr Green Stick Impression Compound 15 Sticks - 113g (00444)



\$65.70 ea + GST

Renfert C&B GEO Classic Modelling Wax Range 1 x 75g Block

Options Available

From
\$25.93 ea + GST



Barrier Sleeves and Film Roll



Syringe Sleeves
500pc

\$16.20 + GST



Blue Film Roll
1200 Sheets

\$22.50 + GST



ADD TO CART

ONLINE ORDER FREIGHT OFFERS

Free Shipping or *Subsidised Flat Rate Shipping Offers* Available
for most orders valued over \$250.00 ex GST
with 2 or more items when placed Online through our Web Shop

*Conditions Apply

Excludes some items (eg. Stone, Equipment, Large Items) and delivery locations.
It will tell you at checkout if promo applies.
Freight promo not available for phone orders.

Valid SEPT / OCT 2024 - All Prices Exclude GST - Prices Do Not Include Freight - Durodent Reserves The Right To Amend Errors and Prices - Subject to stock availability

ONLINE STORE
www.durodent.com.au

SAVE TIME - ORDER ONLINE

ADD TO CART

Speak to our team:
Call us 1800 653 439

DURODENT ACRYLIC TEETH RANGE

durodent
dental



KULZER
MITSUBISHI CHEMICALS GROUP

DELARA 3 Layer Modern Aesthetics

Everybody's Delara – a modern all-round tooth line for everyday use with lifelike aesthetics.

A single tooth line that meets every need of dental laboratories – delivering aesthetics, form and function, with a lively surface structure and unique layering, along with modern and bright shading for a fresh and vivid appearance.

3 Layer Incomp Cad Cam
Shades Available: A1 – D4 + 2 Bleach Shades BL2 / BL3
Moulds: 18 Upper / 5 Lower Anterior Moulds - 8 Posteriors Moulds

Only
\$14.85

Per Card or Less

BUY 50-99
SAVE 10%
\$13.36 ea

BUY 100+
SAVE 15%
\$12.62 ea

INTRO CASE +
CUSTOM BULK
OFFERS
AVAILABLE



Kulzer Pala MONDIAL Acrylic Teeth – CAD CAM Designed Perfection – Wear Resistance Tech

Superior Wear Resistance Properties For Fixed or Removable Prosthetics
Highly Pleasing Classic Aesthetics
NanoPearl Technology – CAD CAM Design makes setup a breeze

Used extensively for All On 4 and Implant Cases

The Pala Mondial tooth line with its wide range of shapes offers a high degree of flexibility in all areas of removable and fixed prosthetics.

Anteriors
RRP \$40.60 ea
Posteriors
RRP \$34.45 ea
Per Card or Less

BUY
30-49
SAVE 10%

BUY 50 +
SAVE 15%

100 + BULK
DISCOUNTS +
INTRO CASE +
CUSTOM BULK
OFFERS
AVAILABLE



Special Intro Cases
Limited Available
A2 / A3 or A3.5 Shades
*New customer only -
Limit 1

Kulzer IDEALIS 8 Posteriors – Flat Cusp

5 Deg Flat Cusp Posteriors
Superior Wear Resistance Properties

Ideal posteriors for implants and prostheses for the elderly patients or specific cases requiring low cusps. CAD CAM Manufactured

Perfectly matched with Kulzer Premium,
Mondial and Delara Anteriors



Posteriors
RRP \$39.41 ea
Per Card or
Less

Kulzer Pala PREMIUM 6 / 8

For The Highest Demands - Mix and Match with Mondial

Pala Premium is the tooth line for the highest demands for aesthetics, function, and natural appearance.
Featuring Kulzer nano pearl technology
Removeable and Fixed Implant Retained Dentures

With their fully anatomical morphology and natural light dynamic, Premiums yield the very highest results, particularly for removable and fixed dentures, and thus facilitates effective integration to the residual dentition.

Shades: 16 VitaPan Shades (A1 – D4), 2 Bleach Shades (BL2-BL3)



Anteriors
RRP \$55.35 ea
Per Card or Less



Posteriors
RRP \$45.83 ea
Per Card or Less

Valid SEPT / OCT 2024 - All Prices Exclude GST – Prices Do Not Include Freight - Durodent Reserves The Right To Amend Errors and Prices – Subject to stock availability

ONLINE STORE
www.durodent.com.au

SAVE TIME - ORDER ONLINE



Speak to our team:
Call us 1800 653 439